



- AlternativeTurorialbox
- Today we will use 1 Eyetracker!
- Create account at sr-support

## Today


- 1. Short: How does the Eyelink 1k+ work?
- 2. How to do gaze control in PTB & Psychopy?
- 3. Preparing EEG/ET coregistration in your experiment
- 4. Storing behavioral data in your eyetrack
- 5. Proper calibration
- 6. Converting .edf files: stupid (Matlab) & smart (R)
- 7. Preprocessing edf files in R

Eyelink



How does the EL1000+ work?






read. The apparatus seemed to work equally well under all the various conditions, even when the speed reached an average of twelve words per second. In order to


- Infrared based eyetracker
- Infrared light is emitted to the eye
- The lens records the eye from another angle
- EL-host-PC detects pupil and Purkinje image

pupil



Purkinje images

#### Calibration



#### calibration


- Calibration creates constellation-dependent coordinate system
- Validation takes similar measurements and checks how much the screen position differs from the system's prediction



Experiments



Matlab / python



# Install this (to record)



- SR Developer's pack
- Psychtoolbox comes with Eyelink toolbox
- Psychopy comes with \*too outdated\* Eyelink module
  - See my readme on github for how to update (don't use pip!)
- My set of high-level Eyelink functions
  - github.com/wanjam/wm\_utilities/Eyelink

# Install this (to analyze)

<u>===</u> พพเ

- SR Developer's pack
- Matlab:
  - Coregister EEG/ET: Olaf Dimigen's <EYE-EEG>
    - Eeglab plugin manager is outdated. Use github: olafdimigen/eye-eeg
- R:
  - Parse edf files: Jason Hubbard's <edfr>
    - devtools::install\_github('jashubbard/edfr')
  - Massage data: Hedderik van Rijn's <pR>, my <wmR>
    - devtools::install\_github('wanjam/wm\_utilities/helper-funs/wmR')
    - devtools::install\_github('hedderik/pR')

## SR-dev kit

- SR research provides a "developer's kit"
- Newest version always in forum
- Ubuntu: PPA (read instructions \*very\* careful)
- Windows/Mac: zip/dmg
- Core of everything (written in C++) → install!
- Have a look: sr-support.com










≡ wwu

Experiments



## Functions for Experiments



#### Readme online! https://github.com/wanjam/wm\_utilities/Eyelink

included here with SR-Research's permission. For the most recent version check their support forum.

#### <sup>™</sup> Usage

Note that this is probably quite specific to our lab. Use this in any other lab at your own risk.

#### Matlab

The Matlab code consists of five simple functions: EyelinkStart.m, EyelinkStop.m, EyelinkRecalibration.m, EyelinkgGetGaze.m, and EyelinkAvoidWrongTriggers.m


A typical workflow in Psychtoolbox would look like this:

1. Start your Experiment with the typical settings until you open a window. Typically, you'll have a line like this: wPtr = Screen('Openwindow', P.PresentScreen, [ 0.5, 0.5, 0.5 ],[],[],[],[],1);



## wm\_ utilities


#### I'll only talk about Matlab – Same functions for Python!




∨ ひ "helper\_funs...

## wm\_ utilities

### I'll only talk about Matlab – Same functions for Python!





r PC > Desktop > git-repos > AlphaFix > helper\_funs

#### Connect

- [P, wPtr] = EyelinkStart(P, wPtr, "filename")
   doc EyelinkStart
- P = your structure with parameters can be empty, can contain extra config stuff

#### Always back assign P!

- Starts and configures connection between EL and PC
- Opens a file to write to
- Defines sampling rate and what to store
- Defines locations of calibration points
- Nulls the TTL port
- Runs an initial calibration



## Recalibrate

### EyelinkRecalibration(P)

Recalibrates the Eyetracker

#### Useful when:

- Break
- Subject moved a lot and eye is lost

## Control Fixation

[didrecal, FixationOnset, hsmvd] =
EyelinkControlFixation(P, Tmin, Tmax, loc,
maxDegDeviation, eyelinkconnected, pixperdeg, dorecal,
Igoreblinks)

Highest possible level in wm\_utilities

- Takes a target location (loc)
- Subject has <Tmax> seconds to look at the target
- Subject has to look at target for <Tmin> seconds
- If not: start a recalibration (If you want that)
- Do you want to ignore blinks?

## Get Gaze

[x, y, hsmvd] = EyelinkGetGaze(P, Ignoreblinks,
OverSamplingBehavior, targetXY,
FixLenDeg, eyelinkconnected, pixperdeg)

Lower level: Just get the gaze position and check whether that differs from the last position.

- Checks if the Eyetracker has a new sample available
- If so, return the (x, y) pixel coordinate pair of the current gaze
- In addition, check if the current gaze differs more than <FixLenDeg> degree from targetXY. If so, <hsmvd> is true (1).

## Behavioral data

#### [message] = EyelinkSendTabMsg(varargin)

stores behavioral data coregistered with time and ttl into the .edf file

- Takes as many input arguments as you want
- Converts everything to strings
- Concatenates everything tab-delimited
- Prepends a "<"
- Sends <message> to .edf file

#### Stop recording


#### EyelinkStop(P, eyelinkconnected)

Finishes recording and puts Eyelink to standby

- Stops the recording
- Copies the final file from the EL-host PC to the stimulation PC
- Then puts the Eyelink to standby

## Status Notifier

EyelinkNotify(message, eyelinkconnected)





#### EyelinkAvoidWrongTriggers

Prevents EEG-trigger-mess
Should be at the top of your <u>non-eyetracking</u> experiment

- Checks if the Eyelink-host pc is turned on
- Throws an error if that's the case



## Parameters

,eyelinkconnected' = boolean – if false, all functions are dummies Dummymode doesn't always work for me, but it's worth trying



Implement Gaze control

- Open a terminal and run >ptb3-matlab
- Go to /home/tutorialbox/Documents/MATLAB/Tutorial/Task
- Start the experiment once by running
  - >TutorialExperiment.m

- Add code to start and stop the Eyelink
- Display "Ruebezahl" as status message, right after EyelinkStart
- Display the current trial number as status message
- Recalibrate after each break.
- Control that subjects look at the fixation symbol.
- Send "triggers" (i.e., messages) after all relevant flips.
- Store the resulting behavioral data in the .edf
- If time: try to make the target's position gaze-dependent



Analyze



## Preprocessing



- Coregistration with EEG data:
  - Fixation Related Potentials → follow EYE-EEG/Elektro-Pipe
- Pupil Dilation
- Saccades
- Fixations
- Gazepaths
- ...

→ Most stuff possible after following this template...

## Parse Matlab

- edf2asc: binary available in the terminal (from the devpack)
- EYE-EEG
  - $parseyelink(file.asc) \leftarrow converts .asc to .mat$
  - pop\_importeyetracker(EEG, file.asc) ← coregistration
- Elektro-Pipe
  - func\_importEye combines all three and plots quality
  - Called by prep01\_preproc

### Parse R

- Jason Hubbard's edfr uses the devpack's C++ code directly
  - No need for edf2asc
  - Most sophisticated package for parsing (afaik)
  - edf.all(), edf.recordings(), edf.events()
  - = ET-data, Header, TTL

- Open Rstudio in the VirtualBox
- You'll find a file called tutorial\_preprocEDF.R
- Let's first walk through that...

- Challenge:
  - Set the output time vector relative to fixation onset
  - Go to line 213. Try to add more info (as added in the experiment)
  - Add the new info to all necessary lines.